

IMPLANTOLOGY

motor system by Bien-Air



PURE SIMPLICITY

S W I S S  M A D E



NEW CHIROPRO

IMPLANTOLOGY
motor system

NEW CHIROPRO

PURE SIMPLICITY

IMPLANTOLOGY
motor system

Simplicity. The cornerstone of the new implant motor developed by Bien-Air Dental.

When designing this new system, everything was geared towards making this ingenious solution simpler to use, down to the smallest detail.

The **New Chiropro** system is not only simple to use, but it is your ideal solution for precise and smooth implant placement.

SWISS  MADE

ABSOLUTE PRECISION

Thanks to its high-precision electronic system, the torque selected on the console corresponds exactly to the torque obtained at the output of your instrument*. No further calibration is needed. It's as simple as that.

PRE-SET SURGICAL PROTOCOLS

Thanks to the pre-set surgical protocols, you can carry out your procedures in no time. Use your patient's clinical data to determine the number of sequences required to perform your procedure - 3, 4 or 5. The parameters of each step are then automatically displayed. You can now simply treat your patient.

CLEAR, CONCISE INTERFACE

The simple and discreet screen provides you with access to all the information you need to ensure your procedure is a success: choice of handpiece, speed, torque, irrigation flow and direction of rotation.

INTUITIVE NAVIGATION

With its single rotary knob, navigation is easier than ever before. Simply turn the knob to scroll through the menus, then click to confirm your choice.

AND MAINTENANCE IS SIMPLE, TOO

The single rotary knob is easy to remove from the unit for sterilization. In addition, the control panel is protected by highly resistant glass that is easy to clean.

* According to tests carried out in Bien-Air Dental's research laboratory, the torque values displayed on Bien-Air Dental systems are considerably more accurate (+/- 5% margin of error) than those displayed on competitors' systems (+/- 10% margin of error).

NEW CHIROPRO

General dental practitioner or experienced surgeon, are you seeking a fast, simple and efficient solution for your implant dentistry procedures?

Look no further than the **Chiropro**. Your perfect solution for precise and smooth placement of implants.

IMPLANT DENTISTRY

Powered by the **Chiropro**, the new MX-i micromotor combined with the CA 20:1 handpiece offer you the very best rotary technology for all your implant-related procedures.

Coupled with the MX-i micromotor, the CA 20:1 handpiece (also available in Kirschner/Meyer and Micro-Series versions) provides unparalleled speed stability for precise and smooth placement of implants. As the most durable implantology handpiece on the market*, the CA 20:1 handpiece has a unique internal irrigation system. Thus, you are not bothered by the irrigation line when handling the instrument. You can carry out your procedures easily and with confidence.



PURE SIMPLICITY

IMPLANTOLOGY
motor system



S W I S S  M A D E

* According to tests carried out in Bien-Air Dental's research laboratory, the drive system of the CA 20:1 CA handpiece has been optimized to withstand more than 1600 uses at 70 Ncm. This substantially increases the instrument's service life.

NEW CHIROPRO

NEW CHIROPRO



NEW CHIROPRO + MAINTENANCE TRIOPACK

1 x SET CHIROPRO 3RD GEN REF 1700708-001



- ✓ Intuitive navigation using a single rotary knob
- ✓ Micromotor MX-i LED with ceramic ball bearings 3-year warranty
- ✓ Maximum torque 70 Ncm

Pack:

- 1 x Chiropro 3rd Gen Console
- 1 x MX-i LED 3rd Gen Micromotor, 3-Year Warranty
- 1 x MX-i LED 3rd Gen Cable
- 1 x Foot Pedal
- 1 x Handpiece holder
- 5 x Irrigation Lines
- 2 x Sterile Protective Films
- 1 x Trio Pack Maintenance

NEW CHIROPRO with 1 Contra-angle CA 20:1 L + MAINTENANCE TRIOPACK

1 x SET CHIROPRO 3RD GEN REF 1700707-001



- ✓ Intuitive navigation using a single rotary knob
- ✓ Contra-angle with light CA 20:1L Micro-Series 2-year warranty
- ✓ Micromotor MX-i LED with ceramic ball bearings 3-year warranty
- ✓ Maximum torque 70 Ncm

Pack:

- 1 x Chiropro 3rd Gen Console
- 1 x Contra-angle CA 20:1 L Micro-Series
- 1 x MX-i LED 3rd Gen Micromotor, 3-Year Warranty
- 1 x MX-i LED 3rd Gen Cable
- 1 x Foot Pedal
- 1 x Handpiece holder
- 5 x Irrigation Lines
- 2 x Sterile Protective Films
- 1 x Trio Pack Maintenance

NEW CHIROPRO with 2 Contra-angles CA 20:1 L + MAINTENANCE TRIOPACK

1 x SET CHIROPRO 3RD GEN REF 1700707-001
+ 1 x CA 20:1 L MICRO-SERIES REF 1600692-001



- ✓ Intuitive navigation using a single rotary knob
- ✓ Contra-angle with light CA 20:1L Micro-Series 2-year warranty
- ✓ Micromotor MX-i LED with ceramic ball bearings 3-year warranty
- ✓ Maximum torque 70 Ncm

Pack:

- 1 x Chiropro 3rd Gen Console
- 2 x Contra-angle CA 20:1 L Micro-Series
- 1 x MX-i LED 3rd Gen Micromotor, 3-Year Warranty
- 1 x MX-i LED 3rd Gen Cable
- 1 x Foot Pedal
- 1 x Handpiece holder
- 5 x Irrigation Lines
- 2 x Sterile Protective Films
- 1 x Trio Pack Maintenance

INSTRUMENTS AND ACCESSORIES

INSTRUMENTS AND ACCESSORIES

**CONTRA-ANGLE
WITH LIGHT CA 20:1 L**
1 x CONTRA-ANGLE CA 20:1 L MICRO-SERIES
REF 1600692-001
OR
1 x CONTRA-ANGLE CA 20:1 L STANDARD
REF 1600598-001



- Patented external irrigation technology
- Multi-strand optical glass technology
- Highly resistant to corrosion

**PM 1:1
EXTERNAL SPRAY**
1 x PM 1:1 REF 1600052-001



- External spray
- Compatible with Micro-Series micromotors
- Without light

**KIT MOTOR MX-i
+ CABLE
(FOR NEW CHIROPRO 3rd GEN)**

SET MX-i LED + CABEL MX-i LED REF 1700744-001
WITH:
• 1 x MICROMOTOR MX-i LED 3rd GEN REF 1601008-001
• 1 x CABLE MX-i LED 3rd GEN REF 1601009-001



**ROTARY KNOB
(FOR NEW CHIROPRO 3rd GEN)**

1 x ROTARY KNOB
REF 1307031-001



Sterilizable rotary knob

**IRRIGATION
LINES**

10 x IRRIGATION LINES
REF 1500984-010
10 x IRRIGATION LINES
KIRSCHNER/MEYER
REF 1501635-010



Set of 10 irrigation lines
sterile single use

S W I S S  M A D E

WWW.BIENAIR.COM

**STERILE
PROTECTIVE FILMS**

10 x STERILE PROTECTIVE FILMS FOR IPAD
REF 1501746-010
or 10 x STERILE PROTECTIVE FILMS
FOR CHIROPRO 3rd GEN
REF 1502329-010



Set of 10 sterile protective films
for single use

MAINTENANCE

MAINTENANCE

LUBRICARE 2 + MAINTENANCE DUOPACK

1 x LUBRICARE 2 REF 1601156-001



Lubricare 2 is compatible with all Bien-Air's turbines and electrical handpieces, Lubriflu-id® and Spraynet®, as well as with hand-pieces and cleaning or lubrication products from other manufacturers. Addi tional couplings and spray holder may be needed to use products not supplied by Bien-Air.

Compatible with handpieces
and cleaning products
from other manufacturers

LUBRIMED
(BOX OF 6 CARTRIDGES)
1 x LUBRIMED REF 1600037-006



DUO PACK
SPRAYNET / LUBRIFLUID
1 x LUBRIFLUID
REF 1600064-001
1 x SPRAYNET
REF 1600036-001



TRIO PACK
SPRAYNET / LUBRIFLUID / AQUACARE
1 x LUBRIFLUID
REF 1600064-001
1 x SPRAYNET
REF 1600036-001
1 x AQUACARE
REF 1600617-001





PURE SIMPLICITY

NEW CHIROPRO

WWW.BIENAIR.COM