

Hillrom is a part of Baxter

Baxter



Hillrom[™]

Welch Allyn Green Series
Medical Exam Lighting



BRIGHT LIGHT WHEREVER YOU NEED IT, FOR LESS

Our **Green Series** Medical Exam Lights give off a brighter, whiter LED light and consume less energy than traditional halogen lights. See true tissue color in a well-defined lit field while also reducing your cost of ownership.

Welch Allyn **Green Series** Medical Exam Lights



Green
Series Exam
Light IV

GREEN SERIES EXAM LIGHT IV

- Ideal for OB/GYN offices or specialties that require focused light in specific areas
- One LED with 3X the typical light output of halogen
- Get pinpoint accuracy with a compact, unobtrusive head and a long, flexible light pipe



GREEN SERIES 300 GENERAL EXAM LIGHT

- Ideal for primary care facilities, emergency rooms, outpatient facilities and dermatology
- One LED gives a perfect combination of intense light output and cool operation, all delivered in a compact design
- Affordable price and rugged performance

Green Series
300 General
Exam Light



NOW AVAILABLE

Midmark accessory mounting bracket for Green Series Exam Light IV, 300, and 600 models.

#48956

BRIGHTER

- View true tissue color with 5,500° K
- Get a focused, uniform light with crisp edges and no dark or hot spots



GREENER

- Lower costs with long-lasting LEDs that do not have to be replaced
- Consume less energy

MORE PATIENT-FRIENDLY

- Reduce risk of cross-contamination with touchless controls and Welch Allyn Disposable Sheaths
- Improve patient comfort and safety with minimal heat dissipation



Welch Allyn Green Series Procedure Lights



Green Series
600 Minor
Procedure Light

GREEN SERIES 600 MINOR PROCEDURE LIGHT

- Ideal for minor procedures and exams in any healthcare environment
- Three LEDS and a wider head give exceptional light in demanding viewing conditions
- Adjust for optimal output with touchless intensity controls

GREEN SERIES 900 PROCEDURE LIGHT

- Ideal for procedures in your hospital, ambulatory care center, outpatient facility or physician's office
- Six LEDs with intensity controls and outstanding durability, maneuverability and spot quality
- Maximize floor space and get direct light wherever you need it with ceiling mount option



Green Series 900
Procedure Light

Welch Allyn Green Series Medical Headlights



GREEN SERIES PORTABLE HEADLIGHTS

- Ideal for general examination, minor suturing and foreign body removal
- Enjoy maximum comfort with a weight of only 4 ounces and choice of rigid or terrycloth headband
- Rechargeable batteries within headband give 60 minutes of continuous runtime from full charge
- LED provides up to 10,000 hours of life
- Coaxial luminaire gives well-defined, shadow-free illumination



GREEN SERIES PROCEDURE HEADLIGHT

- Ideal for any procedure or application that requires a small, bright light, such as foreign body removal, OB/GYN and dermatology
- Enjoy maximum comfort with a weight of only 6 ounces
- Use with direct power supply or a rechargeable battery that gives 60 minutes of continuous runtime from full charge
- LED provides up to 50,000 hours of life
- Coaxial luminaire gives well-defined, shadow-free illumination

When you need targeted, well-defined illumination, look to our high quality, long-lasting and reliable **Green Series** Procedure Headlights to improve staff satisfaction and enhance patient care.



Specifications



EXAM LIGHTS

PROCEDURE LIGHTS

HEADLIGHTS

	Green Series Exam Light IV	Green Series 300 General Exam Light	Green Series 600 Minor Procedure Light	Green Series 900 Procedure Light	Green Series Portable Headlights	Green Series Procedure Headlight
Number of LEDs	1 White LED	1 White LED	3 White LEDs	6 White LEDs	LED	5,700° K
Color Temperature	5,500° K	5,500° K	5,500° K	5,500° K	5,300° K	5,700° K
Light Intensity at Typical Working Distance	20,000 lux (1,858 footcandles)	12" distance: 30,000 lux (2,787 footcandles) 20" distance: 15,000 lux (1,393 footcandles)	High Intensity Setting: 12" distance: 75,000 lux (6,967 footcandles) 20" distance: 40,000 lux (3,716 footcandles) Low Intensity Setting: 12" distance: 55,000 lux (5,109 footcandles) 20" distance: 30,000 lux (2,787 footcandles)	High Intensity Setting: 24" distance: 100,000 lux (9,290 footcandles) 40" distance: 55,000 lux (5,109 footcandles) Low Intensity Setting: 24" distance: 75,000 lux (6,967 footcandles) 40" distance: 40,000 lux (3,716 footcandles)	2,000 lux, 40 lumens	9,000 lux, 120 lumens
LED Life	50,000 hours	50,000 hours	50,000 hours	50,000 hours	10,000 hours	18 in (46 cm)
Working Distance	16 in (40 cm)	12 in (30.5 cm) to 20 in (50 cm)	12 in (30.5 cm) to 20 in (50 cm)	24 in (60.9 cm) to 40 in (100 cm)	16 in (40 cm)	18 in (46 cm)
Spot Size at Working Distance	4 in (10 cm) to 10 in (25 cm)	5 in (13 cm) to 6 in (15 cm)	6 in (15 cm) to 7 in (18 cm)	9 in (23 cm) to 10 in (25 cm)	4 in (10 cm) to 8 in (20 cm)	1 in (2.5 cm) to 6 in (15 cm)
Electrical	Input: 100-240 V~ 50 Hz-60 Hz, 0.30 A max	Input: 100-240 V~ 50 Hz-60 Hz, 0.35 A max	Input: 100-240 V~ 50 Hz-60 Hz, 0.65 A max	Input: 100-240 V~ 50 Hz-60 Hz, 1.00 A max	100-240 V~ 50 Hz-60 Hz, 0.2 A max	100-240 V~ 50 Hz-60 Hz, 0.2 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long		
Weight	Total Table/Wall Mount: 5.30 lbs (2.40 kg) Total Mobile Stand: 12.40 lbs (5.63 kg)	Total Table/Wall Mount: 5.65 lbs (2.56 kg) Total Mobile Stand: 12.75 lbs (5.78 kg)	Total Table/Wall Mount: 6.45 lbs (2.93 kg) Total Mobile Stand: 17.35 lbs (7.87 kg)	Total Wall Mount: 17.90 lbs (8.12 kg) Total Mobile Stand: 44.85 lbs (20.34 kg) Total Ceiling Mount: 4.50 lbs (11.11 kg)	4 ounces (113 g)	6 ounces (170 g)
Mounting Options	Table/Wall Mount Mobile Stand Compatible with Midmark exam lighting tables with Adapter plate	Table/Wall Mount Mobile Stand Compatible with Midmark exam lighting tables with Adapter plate	Table/Wall Mount Mobile Stand Compatible with Midmark exam lighting tables with Adapter plate	Wall Mount Mobile Stand Ceiling Mount	Head Mount–Terrycloth or Rigid Headband	Head Mount–Rigid Headband
Accessory Options	Disposable Sheath #52640	NA	Disposable Sheath #52630	Disposable Sheath #52630	Soft Case #05232-U	Carrying Case #49099, Battery Pack with Charger #75260
Warranty	5-year limited	5-year limited	5-year limited	5-year limited	1-year limited	1-year limited

CHOOSE FROM A VARIETY OF MOUNTING OPTIONS



PRODUCT AND ACCESSORY ORDERING INFORMATION

GREEN SERIES EXAM LIGHT IV

48810: Green Series Exam Light IV with Table/Wall Mount
48950: Mobile Stand for Green Series Exam Light IV
52640: Disposable Sheath for Green Series Exam Light IV (25/box, 5 boxes/case)
48956: Adapter plate for Midmark exam light mounting brackets

GREEN SERIES 300 GENERAL EXAM LIGHT

44400: Green Series 300 General Exam Light with Mobile Stand
44410: Green Series 300 General Exam Light with Table/Wall Mount
44215: 12" Wall Mount Extension
48960: Heavy Duty/Tall Mobile Stand
48956: Adapter plate for Midmark exam light mounting brackets

GREEN SERIES 600 MINOR PROCEDURE LIGHT

44600: Green Series 600 Minor Procedure Light with Mobile Stand
44610: Green Series 600 Minor Procedure Light with Table/Wall Mount
44215: 12" Wall Mount Extension
48960: Heavy Duty/Tall Mobile Stand
52630: Disposable Sheath for Green Series 600/
Green Series 900 (50/box, 5 boxes/case)
48956: Adapter plate for Midmark exam light mounting brackets

GREEN SERIES 900 PROCEDURE LIGHT

44900: Green Series 900 Procedure Light with Mobile Stand
44900-C: Green Series 900 Procedure Light with Ceiling Mount
44900-W: Green Series 900 Procedure Light with Wall Mount
52630: Disposable Sheath for Green Series 600/
Green Series 900 (50/box, 5 boxes/case)

GREEN SERIES PORTABLE HEADLIGHT

46070R: Green Series Portable Headlight with Rigid Headband and Direct Power Supply/Charger
46070: Green Series Portable Headlight with Terrycloth Headband and Direct Power Supply/Charger

GREEN SERIES PROCEDURE LIGHT

49020: Green Series Procedure Headlight with Direct Power Supply
75260: Green Series Procedure Headlight Battery Pack with Charger

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

Baxter, Hillrom, and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

APR129104 Rev 2 21-NOV-2022 ENG - US